COCID-19(2019-nCOV) Coronavirus Ig G/IgM Rapid - Test kit

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COVID-19 (2019-NCOV) CORONAVIRUS RAPID TEST KIT

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Shenzhen Lvshiyuan Biotechnology Co., Ltd
Add: 101,201,301, D Building, No. 2 Industrial Avenue, Buxin
Village, Buxin Community, Dapeng Subdistrict Office,
Dapeng New District, Shenzhen, 518120 China

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COVID-19 (2019-nCoV) Coronavirus IgG/IgM Rapid Test Kit
(Colloidal gold)

Shenzhen Lvshiyuan Biotechnology Co.,Ltd.

SHOT ON MI 9
BY PAVVCIO

COVID-19 (2019-nCoV) Coronavirus IgG/IgM Rapid Test Kit
(Colloidal gold)

IVD

SHOT ON MI 9
BY PAVVCIO

Batch NO.: KB200303

MFG : 20200318

Expiry Date: 20210317

Store at :2-30°C

Shenzhen Lvshiyuan Biotechnology Co., Ltd.
Add: 101,201,301, D Building, No. 2 Industrial Avenue, Buxin
Village, Buxin Community, Dapeng Subdistrict Office,
Dapeng New District, Shenzhen, 518120, China.

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
R&D Joint by Third Affiliated Hospital of Sun Yat-sen University
School of Biology and Medicine, Guangdong University of Technology
Shenzhen Lvshiyuan Biotechnology Co.,Ltd

Shenzhen Lvshiyuan Biotechnology Co.,Ltd

101,201,301, D Building, No.2 Industrial Avenue, Buxin Village, Buxin Community,
Dapeng Subdistrict Office, Dapeng New District, Shenzhen, 518120 China

EC REP

Obelis s.a.
Bd Général Wahis 53 1030 Brussels Belgium

Tel: 86-755-28438788
Email: info@lsybt.com

Fax: 86-755-28938800
www.lsybt.com

Declaration of Conformity

According to annex III of the Council Directive 98/79/EC on in vitro diagnostic medical device
We,

NAME: Shenzhen Lvshiyuan Biotechnology Co., Ltd

ADDRESS : 101, 201, 301, D Building, No.2 Industrial Avenue, Buxin Village, Buxin Community, Dapeng Subdistrict Office, Dapeng New District, Shenzhen, 518120 china

Declare under our sole responsibility that the following in vitro diagnostic medical devices other than those covered by annex II and devices for performance evaluation

list of products

COVID-19(2019-nCoV) Coronavirus IgG/IgM Rapid Test Kit

Meet the provisions of the Council Directive 98/79/EC concerning medical devices which apply to them.

Undersigned declares to fulfill the obligations imposed by Annex III section 2 to 5:

- availability of the technical documentation set in Annex III (section 3), allowing the assessment of conformity of the product with the requirements of the Directive.
- the manufacturer shall take necessary measures to ensure that the manufacturing process follows the principles of quality assurance as appropriate for the products manufactured (Annex III section 4).
- the manufacturer shall institute and keep up to date a systematic procedure to review experience gained from devices in the post-production phase and to implement appropriate means to apply any necessary corrective actions (Annex III section 5).

Conformity assessment was performed according to Article 9 (7) and Annex III, section 3.

Our current Quality System is formatted to international standards:

**PLEASE SPECIFY
ISO9001:2015**

Corporate Contact Information

COMPANY NAME: Shenzhen Lvshiyuan Biotechnology Co., Ltd

COMPANY ADDRESS 1: 101, 201, 301, D Building, No.2 Industrial Avenue, Buxin Village, Buxin Community, Dapeng Subdistrict Office, Dapeng New District, Shenzhen

COMPANY ADDRESS 2 : NA

COMPANY PHONE: 18922879555

COMPANY FAX: NA

COMPANY EMAIL: szlsy2003@163.com

RESPONSIBLE PERSON'S name: Jiang Yongqing

Position: Vice General manager

SIGNATURE :

Date :

Stamp

European Authorized Representative:

Registered Address:

Obelis s.a.

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B-1030 Brussels, Belgium

Phone: 32.2.732.59.54

Fax: 32.2.732.60.03

E-mail: mail@obelis.net

Representative: Mr. Gideon ELKAYAM (CEO)

CERTIFICATE OF REGISTRATION

The Quality Management Systems of

Shenzhen Lvshiyuan Biotechnology Co., Ltd

Unified Social Credit Code: 914403007576264357

Registration address: 101, 201, 301, D Building, No.2 Industrial Avenue, Buxin Village, Buxin Community, Dapeng Subdistrict Office, Dapeng New District, Shenzhen
Production address: D Building, National Biological Industrial Park Of Marinelife, Binhai No.2 Road, Dapeng, Shenzhen

has been assessed by GIC and complying with

GB/T19001-2016/ISO9001:2015

For the following activities

**Research and development, production and service of
food safety testing kits, animal disease diagnostic kits
and test cards**

Date of Issue: 13 February 2019

Date of Expiry: 12 February 2022

Date of Initial Certification: 13 February 2019

Certificate No.: J19Q2GZ8012523R0M

Scan for certificate status

The granting of this certificate does not mean that the certificate holder can avoid any legal obligation. If the products or activities covered in the scope of certification require administrative license, the certificate shall be only valid within the scope of administrative licensing. The registered organization shall be subject to regular annual supervision by GIC, and the continual validity of the certificate is base upon conformity of audit. Please scan two-dimension code at left to find the certificate information. This certificate can be queried at Certification and Accreditation Administration of the People's Republic of China official website (www.cnca.gov.cn) & GIC website (www.gicg.com.cn)

GIC WeChat public number

Signature: *Helen*

Guardian Independent Certification Ltd

Registered in England
Sovereign House 212-224 Shaftesbury Avenue London England WC2H 8HQ
Accredited by Member of IAF MLA
JAS-ANZ registration no. 335400064, www.jas-anz.org/auk

CERTIFICATE OF IVD NOTIFICATION

Ref. No.: BS 8841-2020

BELGIUM

Order No.: BS 8825-2020

Date: 20/03/2020

THIS IS TO CERTIFY THAT, ACCORDING TO THE COUNCIL DIRECTIVE 98/79/EC, OBELIS S.A. (O.E.A.R.C.) PERFORMED ALL NOTIFICATION DUTIES AND RESPONSIBILITIES AS THE EUROPEAN AUTHORIZED REPRESENTATIVE (EC REP) OF:

NAME: SHENZHEN LVSHIYUAN BIOTECHNOLOGY CO., LTD.

ADDRESS: 101, 201, 301, D BUILDING, NO. 2 INDUSTRIAL AVENUE, BUXIN VILLAGE, BUXIN COMMUNITY, DAPENG SUBDISTRICT OFFICE, DAPENG NEW DISTRICT, SHENZHEN, 518120, CHINA

AS STIPULATED AND DEMANDED BY THE AFOREMENTIONED DIRECTIVE.

The Manufacturer declares that the IVD devices comply with the Directive including all essential requirements.

The Manufacturer has provided Obelis s.a. (O.E.A.R.C.) with all the appropriate declarations according to the 98/79/EC Directive – article 10 requirements including the EC Declaration of Conformity confirming that his In-Vitro Diagnostics medical devices, as stipulated here above, are fulfilling the applicable requirements of the European Council Directive 98/79/EC

The notification of the following In-Vitro Diagnostic medical devices has been completed by Obelis s.a. (O.E.A.R.C.) on the 19/03/2020 in compliance with the European Council Directive 98/79/EC - article 10 requirements.

IN-VITRO DIAGNOSTIC MEDICAL DEVICES: PLEASE SEE ANNEX A - LIST OF DEVICES (1 PAGE: 1 DEVICE)

As of the 20/03/2020, and as long as the manufacturer will continue complying with the hereabove mentioned requirements* he therefore:

- Is required to affix the CE marking on these devices;
- Place these devices in the Territory of Belgium and/or the other EEA Member States (excluding territories not in alignment with Decision 2010/227/EU).

Mr. G. Elkayam CEO
Obelis sa

Obelis European Authorized Representative Center is a member of the European Association of Authorized Representatives (E.A.A.R.), ISO 9001 : 2015 and ISO 13485 : 2016 certified in accordance to the profession of a European Authorized Representative.

** This Certificate will be automatically void if the notification is rejected by the EU Authorities or upon termination of the EAAR agreement.

Registered Address: Bd. General Wafis 53-1030 Brussels | Registered Office Address: Bd Brand Whitlock 30, B-1200 Brussels - Belgium
T. + 32 (0) 2 732 5954 | F. + 32 (0) 2 732 8005 | Email: rsa@obelis.net | Website: www.obelis.net
V3 - ID : 00454718 - 22/02/2019

1

医疗器械产品出口销售证明
CERTIFICATE FOR EXPORTATION OF MEDICAL
PRODUCTS

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Product(s): COVID-19(2019-nCoV) Coronavirus IgG/IgM Rapid Test Kit

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Address of manufacturer: 101, 201, 301, D Building, No.2 Industrial Avenue,
Buxin Village, Buxin Community, Dapeng Subdistrict Office, Dapeng New
District, Shenzhen

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Manufacturing License(s): 粤食药监械生产许 20142513 号

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This is to certify that the above product(s) are not registered in China. The
exportation of the product(s) is not restricted.

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(This certificate is valid for two years from the date of issuance)

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Manufacturer: Shenzhen Lvshiyuan Biotechnology Co., Ltd