



**RANVOO PH3 TEENAGER |  
KOŃCÓWKA DO SZCZOTECZKI  
ZAWIERAJĄCA FLUOR |  
WŁOSIE MIĘKKIE**



link to the product:

<https://www.batna24.com/gb/p/brush-heads/ranvoo-ph3-teenager-koncowka-do-szczoteczki-zawierajaca-fluor-wlosie-miekkie>

Price: **109.99 zł** gross / 89.42 zł net



Manufacturer: RANVOO

Referention number: DNXGKOJPLRMNHM

**Information:**

The RANVOO PH3 Soft sonic brush head with fluoride provides gentle yet effective cleaning of teeth and interdental spaces. The soft, dense bristles work with sonic micro-bubble technology, helping to remove plaque and food debris and deliver fluoride ions during brushing. It is a good solution for people with sensitive gums, supporting enamel mineralisation and protection against tooth decay. The compact, ergonomic design makes it easier to reach hard-to-access areas. Replacing the brush head every 2-3 months helps maintain hygiene and cleaning effectiveness.

**Full product description**


**Care and enamel protection**

The RANVOO PH3 Soft Sonic Head with fluoride has been designed for gentle yet effective oral cleaning. Soft, dense bristles work with the sonic micro-bubble technology of the RANVOO toothbrush, removing plaque and food debris while delivering fluoride ions during brushing.



### For sensitive gums

The Soft fluoride head is suitable for people with sensitive gums. Fluoride released from the bristles supports enamel mineralisation, helping to reduce the risk of tooth decay and maintain a healthy tooth surface during daily brushing.




### Precise and compact design

The PH3 Soft head has a smaller, ergonomic shape, making it easier to reach hard-to-access areas such as the interdental spaces and the gum line. The soft bristles gently surround the tooth surfaces, while the compact design allows thorough cleaning.



### **Optimal compatibility with the PH3 toothbrush**

The head works with the sonic movements of the PH3 toothbrush, combining a high vibration frequency with the action of micro-bubbles, which enhances the removal of deposits and plaque. Bristles of different lengths increase cleaning effectiveness while ensuring comfortable use.



### **Regular replacement and hygiene**

To maintain performance and hygiene, it is recommended to replace the head every 2-3 months, or sooner if the bristles become deformed or worn. A fresh head supports effective cleaning and maintaining fluoride levels in the mouth.

