



Extralink Beauty Ultrasonic Skin Cleaner 3in1 | Ultrasonic Skin Cleaner |



link to the product:

<https://www.batna24.com/gb/p/body-and-face-care/extralink-beauty-ultrasonic-skin-cleaner-3in1-ultrasonic-skin-cleaner->

Price: **89.90 zł** gross / 73.09 zł net

Manufacturer: EXTRALINK



Referention number: DNXXGMMJKRMMKO

Information:

Extralink Beauty Ultrasonic Skin Cleaner 3in1 is a multifunctional device for at-home facial care, combining cavitation peeling, sonophoresis and EMS microcurrents. 28 kHz ultrasound helps remove dead skin cells, excess sebum and impurities, supporting skin cleansing and smoothing. The sonophoresis mode makes it easier to apply active ingredients from serums or skincare products, while EMS provides a gentle massage that supports skin firmness and facial contours. Clear markings and LED indicators make it easy to use, and the 300 mAh battery and compact design ensure convenience at home and when travelling.

Full product description

3-in-1 facial care

Extralink Beauty Ultrasonic Skin Cleaner 3in1 combines **cavitation peeling, sonophoresis and EMS microcurrents**. Three operating modes allow the device to be adapted to the successive stages of home care: cleansing, applying cosmetics and skin massage.



Cavitation peeling

The device uses **ultrasound at a frequency of 28 kHz**. When used on damp skin, the metal spatula helps remove dead skin cells, excess sebum and surface impurities. The treatment supports skin cleansing and helps smooth the skin.



Support for cosmetic products

The **sonophoresis** function uses ultrasonic vibrations to support the application of active ingredients contained in a suitably selected serum or other skincare product. This helps combine device-based care with the skin hydration and nourishment stage.



EMS microcurrents

The **EMS** mode uses gentle electrical impulses during facial massage. The function supports care aimed at improving skin firmness and emphasising facial contours. The device can be moved over selected areas in accordance with the direction indicated in the user manual.



Simple home routine

Clear mode markings and **LED indicators** make the device easy to use. Peeling should be performed on cleansed and moistened skin, while during sonophoresis a cosmetic product intended for this type of care can be applied. After the treatment, the metal tip should be thoroughly cleaned.



Rechargeable design

The rechargeable power supply and **300 mAh capacity** allow the device to be used cord-free during treatment. The slim housing and cover for the metal spatula make it easy to store the beauty accessory in the bathroom or take it on a trip.

