



EXTRALINK MINI MASSAGE STICK



link to the product:

<https://www.batna24.com/gb/p/massage-pillows/extralink-mini-massage-stick>

Price: **16.46 zł** gross / 13.38 zł net

Manufacturer: EXTRALINK



Referention number: DNXGGMKOGRMJJO

Information:

Extralink Mini Massage Stick is a compact muscle electrostimulator using TENS and EMS technologies, which helps relieve pain, reduce tension and support recovery after work, training or prolonged sitting. It offers 8 operating modes and 19 intensity levels, making it easy to adjust to your own needs. It is suitable for massaging the neck, back, waist, hips, legs and abdomen. It has a clear display, simple operation and a lightweight, portable design. The efficient battery allows for long use, and fast USB charging takes only 30 minutes. The durable, easy-to-clean gel pads can be used repeatedly.

Full product description

Extralink Mini Massager

Extralink Mini Massage Stick is a compact **muscle electrostimulator** using **TENS and EMS** technologies to relieve tension and support the body's recovery. The device helps reduce pain in the neck, back and limbs, offering a professional standard of therapy at home. Thanks to its small size, the **portable massager** is ideal for use while resting, working at a desk or after training.



Key product benefits

- **TENS and EMS technology** - effectively helps relax tense muscle structures and ease discomfort.
- **Versatile action** - the design allows massage of the **neck, waist, hips, legs and abdomen**.
- **Customisable treatment** - choose from **8 operating modes** and as many as **19 intensity levels** of stimulation.
- **Efficient battery** - the **120mAh** battery allows use for up to **32 days** on a single charge.
- **Fast recharging** - a full charging cycle via the **USB** port takes just **30 minutes**.
- **Convenience and hygiene** - durable gel pads can be used up to **60 times** and are easy to clean.



TENS and EMS therapy

The device uses low-frequency electrical impulses in the range of **1-100Hz**, which contribute to the direct stimulation of sensory neurons and muscle fibres. **Electrostimulation** helps block pain signals and stimulates blood circulation, speeding up recovery processes. The **Extralink Mini Massage Stick** is a good choice for people dealing with a sedentary lifestyle or excessive physical exertion, helping to maintain comfort every day.



Intuitive control

Using the massager has been simplified to a minimum, making it accessible to every user. A clear **digital display** shows the currently selected programme, and the control buttons allow smooth adjustment of the settings. Simply attach the **electrostimulator** to clean skin and select the appropriate mode with the "M" button. Power adjustment using plus and minus allows precise matching of the intensity to the user's sensitivity threshold.




Mobility and durability

The housing, made from a combination of **ABS, silicone and PU** materials, guarantees durability while keeping the weight low. The compact **150x70 mm** dimensions mean the massager fits in hand luggage or a handbag. The device offers freedom of movement, as it does not require a constant mains connection while in use. With a voltage of **3.7V** and an optimal output of **35mA**, the massage session is safe and effective wherever you are.



Effective therapy plan

For optimal health results, it is recommended to use the device for **10 minutes a day** in regular cycles lasting at least a week. Consistency makes the muscles more flexible and chronic tension subsides. The **Extralink Massager** automatically ensures session safety, operating only after the pad has been correctly attached to the body. It is a solution that helps maintain physical fitness without the need to visit a physiotherapist.



Economy and ecology

The use of a **rechargeable battery** eliminates the need to constantly buy disposable batteries, resulting in real savings. A single, short 30-minute charge via the standard **USB port** is enough for over a month of daily treatments. The massage pads are designed for repeated contact with the skin - once they lose their adhesion, simply rinse them with water to restore their properties, which significantly extends the life of the entire set.



Specification

Parameter	Details
Material	ABS + silicone + PU
Dimensions	150 x 70 mm
Massage modes	8 modes and 19 adjustable settings
Frequency	1-100 Hz
Battery capacity	120 mAh
Charging time	30 minutes
Operating time	Up to 32 days (at 15 min per day)
Charging connector	USB
Operating voltage	3.7 V
Maximum power	35 mA